

Caerin 1.1 TFA

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage:

Keep away from moisture

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Caerin 1.1 TFA, a host defense peptide derived from the Australian tree frog Litoria's glandular secretions, exhibits anti-proliferative effects on HeLa cells and induces apoptosis. Additionally, it has the potential to augment adaptive T cell immune responses [1].
Targets(IC50)	Apoptosis

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481