

Bufospirostenin A

Chemical Properties


CAS No. : 2097252-84-9

Formula: C₂₇H₄₀O₄

Molecular Weight: 428.6

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	Bufospirostenin A represents the inaugural (spirostanol) compound isolated from an animal source.
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3332 mL	11.6659 mL	23.3318 mL
5 mM	0.4666 mL	2.3332 mL	4.6664 mL
10 mM	0.2333 mL	1.1666 mL	2.3332 mL
50 mM	0.0467 mL	0.2333 mL	0.4666 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481