

Brain Natriuretic Peptide-45, rat TFA

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from moisture
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|--|
| Description | Brain Natriuretic Peptide-45, rat TFA (BNP-45, rat TFA) is an isoform of the rat brain natriuretic peptide extracted from the rat heart, exhibiting significant hypotensive and natriuretic effects [1]. |
| Targets(IC50) | Natriuretic peptide |
| In vivo | Brain Natriuretic Peptide-45, rat (BNP-45, 0.1-2.0 nmol/kg, i.v.), demonstrates significant natriuretic and hypotensive effects in anesthetized spontaneously hypertensive rats (SHR) and Wistar-Kyoto rats (WKY). Higher BNP-45 concentrations are more effective in reducing blood pressure in SHR, while WKY rats display increased sensitivity to its diuretic and natriuretic actions and greater urinary cGMP excretion. Elevated BNP-45 doses result in prolonged blood pressure reduction and urinary cGMP excretion in WKY [1]. |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481