

Bombinin-BO1

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Bombinin-BO1, an antimicrobial peptide from the skin secretions of the Bombina orientalis toad, exhibits activity against E. coli, S. aureus, and Candida albicans with minimum inhibitory concentration (MIC) values of 64, 64, and 128 mg/L, respectively [1].

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481