

Bis-(-)-8-demethylmaritidine

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Bis-(-)-8-demethylmaritidine, a natural alkaloid, is a potent inhibitor of acetylcholinesterase (AChE) and is utilized in Alzheimer's disease research [1].
Targets(IC50)	Others,Cholinesterase (ChE)

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481