

Biotin-PEG-OH (MW 20000)

**Chemical Properties**

CAS No. :

Formula:

Molecular Weight:

Storage: Keep away from direct sunlight  
Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
Actual storage temperature shall be subject to the COA.

**Biological Description**

Description	Biotin-PEG-OH (MW 20000) is a polyethylene glycol (PEG)-based PROTAC linker used for synthesizing PROTACs [1].
Targets(IC50)	Others

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481