

Biotin-PEG-Biotin (MW 20000)

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: **Keep away from direct sunlight**
Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

Description	Biotin-PEG-Biotin (MW 20000) is a bifunctional PEG derivative with biotin labels at each terminal end. Biotin facilitates protein labeling as an enzyme co-factor, while the hydrophilic and low-toxicity nature of PEG imparts water solubility to the compound.
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481