

Biotin-H10

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Biotin-H10 is a potent inhibitor of anterior gradient homolog 2 (AGR2) with a dissociation constant (K _D) of 6.4 nM, effectively reducing the viability of cancer cells [1].
Targets(IC50)	MDM-2/p53

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481