

Binetrakin

Chemical Properties

CAS No. : 207137-56-2

Formula:

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description

Binetrakin (Interleukin 4 (human)) is a pleiotropic cytokine from T lymphocytes that affects multiple cell types like B cells and T cells, and promotes the differentiation of naive helper T cells (Th0 cells) into Th2 cells, with applications in various biological research contexts [1] [2].

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481