

BigLEN(mouse) TFA

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	BigLEN(mouse) TFA is a proSAAS-derived neuropeptide and GPR171 agonist that both regulates food intake in mice and inhibits the release of glutamate onto parvocellular neurons of the paraventricular nucleus, a mechanism contingent on postsynaptic G protein activation.
Targets(IC50)	GPCR

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481