

Bempikibart

Chemical Properties

CAS No. : 2622254-57-1

Formula:

Molecular Weight: 145.75 kDa

Store at low temperature

Storage: Store at -20°C

Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|---|
| Description | Bempikibart (ADX-914) is a fully humanised monoclonal antibody targeting IL-7R α , capable of inhibiting IL-17 and TSLP (Thymic stromal lymphopoietin), and may be employed in studies of cutaneous inflammation and alopecia. |
| Targets(IC50) | IL Receptor, Immunology/Inflammation related |

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|-----------|-----------|
| 1 mM | 0.0069 mL | 0.0343 mL | 0.0686 mL |
| 5 mM | 0.0014 mL | 0.0069 mL | 0.0137 mL |
| 10 mM | 0.0007 mL | 0.0034 mL | 0.0069 mL |
| 50 mM | 0.0001 mL | 0.0007 mL | 0.0014 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481