

## Belrestotug

## Chemical Properties

CAS No. : 2574438-65-4

Formula:

Molecular Weight:

Storage: -20°C for 1 year  
Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Belrestotug (EOS-448), an antagonistic anti-TIGIT (VSIG9, VSTM3) antibody of the human immunoglobulin G1 kappa (hIgG1 kappa) class, exhibits antineoplastic activity [1].
Targets(IC50)	Immunology/Inflammation related

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481