


BAX-IN-1

Chemical Properties

CAS No. : 331261-91-7
 Formula: C₁₆H₁₄N₆O
 Molecular Weight: 306.32
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	BAX-IN-1 is a potential selective inhibitor of Bcl-2-associated X protein (BAX).
Targets(IC50)	Bcl-2 Family,Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.2646 mL	16.3228 mL	32.6456 mL
5 mM	0.6529 mL	3.2646 mL	6.5291 mL
10 mM	0.3265 mL	1.6323 mL	3.2646 mL
50 mM	0.0653 mL	0.3265 mL	0.6529 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481