


## Barzuxetan

## Chemical Properties

CAS No. :	157380-45-5
Formula:	C <sub>26</sub> H <sub>34</sub> N <sub>4</sub> O <sub>10</sub> S
Molecular Weight:	594.63
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Barzuxetan is a potential research tool in oncology [1].
Targets(IC50)	Drug-Linker Conjugates for ADC

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6817 mL	8.4086 mL	16.8172 mL
5 mM	0.3363 mL	1.6817 mL	3.3634 mL
10 mM	0.1682 mL	0.8409 mL	1.6817 mL
50 mM	0.0336 mL	0.1682 mL	0.3363 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481