

Bamadutide

Chemical Properties

CAS No. : 1596343-09-7

Formula:

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Bamadutide (SAR425899) is a potent dual agonist for the glucagon-like peptide-1 receptor and glucagon receptor, enhancing β -cell function and reducing glucose absorption in vivo, thereby improving postprandial glucose control. This compound has potential applications in type 2 diabetes research [1].
Targets(IC50)	Glucagon Receptor

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481