

## Aurora A inhibitor 3

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Aurora A Inhibitor 3 (Compound 5h) effectively inhibits Aurora-A kinase, exhibiting an IC50 of 0.78 $\mu$ M, and demonstrates cytotoxicity against MCF-7 and MDA-MB-231 cell lines with GI50 values of 0.12 $\mu$ M and 0.63 $\mu$ M, respectively [1].
Targets(IC50)	Aurora Kinase

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481