


## ASGPR ligand-1

## Chemical Properties

CAS No. :	1426160-58-8
Formula:	C <sub>61</sub> H <sub>110</sub> N <sub>10</sub> O <sub>27</sub>
Molecular Weight:	1415.58
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	ASGPR ligand-1 is a compound that targets the asialoglycoprotein receptor (ASGPR) and is utilized in various disease research contexts [1].
Targets(IC50)	Others,LYTACs

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.7064 mL	3.5321 mL	7.0642 mL
5 mM	0.1413 mL	0.7064 mL	1.4128 mL
10 mM	0.0706 mL	0.3532 mL	0.7064 mL
50 mM	0.0141 mL	0.0706 mL	0.1413 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481