

Anti-inflammatory agent 59

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Anti-inflammatory agent 58 is characterized by its ability to inhibit IL-1 β with an IC50 value of 2.28 μ M. Furthermore, it reduces pro-inflammatory gene expression, protein secretion, and NF- κ B phosphorylation [1].
Targets(IC50)	Others,Interleukin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481