

## Antibacterial agent 163

## Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

## Biological Description

Description	Antibacterial agent 163 (compound 1), a hydroxyquinoline derivative, potently inhibits methicillin-resistant Staphylococcus aureus (MRSA) [1].
Targets(IC50)	Others,Antibacterial

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481