

## Angiopep-Bim BH3 hydrochloride

### Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Keep away from moisture

Storage:

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description

Angiopep-Bim BH3 hydrochloride, a peptode that penetrates the blood-brain barrier (BBB), may serve as a tool to study the permeability of central nervous system (CNS) therapeutics [1].

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481