

Albinterferon alfa-2B

Chemical Properties

CAS No. : 472960-22-8
Formula:
Molecular Weight:
Storage: -20°C for 1 year
Actual storage temperature shall be subject to the COA.

Biological Description

| | |
|---------------|---|
| Description | Albinterferon alfa-2B (alb-IFN), a recombinant protein, comprises interferon alfa-2b conjugated to human albumin and is utilized in the study of chronic hepatitis C [1]. |
| Targets(IC50) | Others |

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481