


## Agistatin E

## Chemical Properties

CAS No. :	144096-48-0
Formula:	C <sub>11</sub> H <sub>16</sub> O <sub>5</sub>
Molecular Weight:	228.24
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Agistatin E, a pyranacetal derived from Fusarium species, impedes cholesterol biosynthesis [1].
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.3814 mL	21.9068 mL	43.8135 mL
5 mM	0.8763 mL	4.3814 mL	8.7627 mL
10 mM	0.4381 mL	2.1907 mL	4.3814 mL
50 mM	0.0876 mL	0.4381 mL	0.8763 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481