

Adrenomedullin (rat)

Chemical Properties

CAS No. : 161383-47-7

Formula:

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Adrenomedullin (rat), a potent vasodilator peptide, is secreted by both endothelial cells (EC) and vascular smooth muscle cells (VSMC) [1].
Targets(IC50)	Interleukin

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481