


## Acetyltetrapeptide 11

## Chemical Properties

CAS No. :	928006-88-6
Formula:	C <sub>27</sub> H <sub>38</sub> N <sub>4</sub> O <sub>7</sub>
Molecular Weight:	530.61
Storage:	Keep away from moisture Powder: -20°C for 3 years   In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



## Biological Description

Description	Acetyltetrapeptide-11, a bioactive peptide with anti-aging properties, is reported to be utilized as an ingredient in cosmetics [1].
Targets(IC50)	Procollagen C Proteinase

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.8846 mL	9.4231 mL	18.8462 mL
5 mM	0.3769 mL	1.8846 mL	3.7692 mL
10 mM	0.1885 mL	0.9423 mL	1.8846 mL
50 mM	0.0377 mL	0.1885 mL	0.3769 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481