


Acetyl-(D-Phe2,Lys15,Arg16,Leu27)-VIP (1-7)-GRF (8-27)

Chemical Properties

CAS No. :	202463-00-1
Formula:	C150H246N44O38
Molecular Weight:	3273.83
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	Acetyl-(D-Phe2,Lys15,Arg16,Leu27)-VIP(1-7)-GRF(8-27) is a selective antagonist for the vasoactive intestinal peptide 1 (VIP 1) receptor.
Targets(IC50)	PACAP

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.3055 mL	1.5273 mL	3.0545 mL
5 mM	0.0611 mL	0.3055 mL	0.6109 mL
10 mM	0.0305 mL	0.1527 mL	0.3055 mL
50 mM	0.0061 mL	0.0305 mL	0.0611 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481