

AADvac 1

Chemical Properties

CAS No. : 2307465-24-1

Formula:

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	AADvac1 is a biologically active peptide functioning as an Alzheimer's Disease Peptide Vaccine, composed of Tau(294-305) conjugated to Keyhole Limpet Hemocyanin (KLH).
Targets(IC50)	Others

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481