


## A1AT modulator 1

## Chemical Properties

CAS No. :	2406205-61-4
Formula:	C <sub>21</sub> H <sub>23</sub> FN <sub>2</sub> O <sub>3</sub>
Molecular Weight:	370.42
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	A1AT Modulator 1 potently inhibits Z α1-antitrypsin (A1AT) polymerization, exhibiting a pIC <sub>50</sub> of 8.3 [1].
Targets(IC <sub>50</sub> )	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6996 mL	13.4982 mL	26.9964 mL
5 mM	0.5399 mL	2.6996 mL	5.3993 mL
10 mM	0.270 mL	1.3498 mL	2.6996 mL
50 mM	0.054 mL	0.270 mL	0.5399 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E\_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481