


7(R)-7,8-Dihydrosinomenine

Chemical Properties

CAS No. : 65494-41-9
 Formula: C₁₉H₂₅NO₄
 Molecular Weight: 331.41
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
 Actual storage temperature shall be subject to the COA.



Biological Description

Description	7(R)-7,8-Dihydrosinomenine, an isoquinoline alkaloid extractable from <i>Sinomenium acutum</i> (Menispermaceae), exhibits anti-osteoclastogenic properties [1].
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0174 mL	15.0871 mL	30.1741 mL
5 mM	0.6035 mL	3.0174 mL	6.0348 mL
10 mM	0.3017 mL	1.5087 mL	3.0174 mL
50 mM	0.0603 mL	0.3017 mL	0.6035 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481