

5-O-Primeverosylapigenin

Chemical Properties


CAS No. : 2382681-38-9

Formula: C₂₆H₂₈O₁₄

Molecular Weight: 564.49

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	5-O-Primeverosylapigenin, a member of the primeverosylapigenin class, can be isolated from species within the Gentianaceae family [1].
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.7715 mL	8.8576 mL	17.7151 mL
5 mM	0.3543 mL	1.7715 mL	3.543 mL
10 mM	0.1772 mL	0.8858 mL	1.7715 mL
50 mM	0.0354 mL	0.1772 mL	0.3543 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481