


5-Azacytidine 5'-triphosphate

Chemical Properties

CAS No. :	2226-74-6
Formula:	C ₈ H ₁₅ N ₄ O ₁₄ P ₃
Molecular Weight:	484.14
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	5-Azacytidine 5'-triphosphate (5-aza-CMP), a cytidine analog, selectively inhibits the incorporation of [³ H]CTP into RNA during the DNA-dependent RNA polymerase reaction without affecting the incorporation of [³ H]UTP [1].
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.0655 mL	10.3276 mL	20.6552 mL
5 mM	0.4131 mL	2.0655 mL	4.131 mL
10 mM	0.2066 mL	1.0328 mL	2.0655 mL
50 mM	0.0413 mL	0.2066 mL	0.4131 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481