


4-(Benzyloxy)-2-nitrophenol

Chemical Properties

CAS No. :	96315-18-3
Formula:	C13H11NO4
Molecular Weight:	245.23
Storage:	Keep away from direct sunlight Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	4-(Benzyloxy)-2-nitrophenol serves as a synthetic precursor for constructing fluorescent ion indicators [1].
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	4.0778 mL	20.389 mL	40.778 mL
5 mM	0.8156 mL	4.0778 mL	8.1556 mL
10 mM	0.4078 mL	2.0389 mL	4.0778 mL
50 mM	0.0816 mL	0.4078 mL	0.8156 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481