


## 3-Hydroxyirisquinone

## Chemical Properties

CAS No. :	100205-29-6
Formula:	C <sub>24</sub> H <sub>38</sub> O <sub>4</sub>
Molecular Weight:	390.56
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	3-Hydroxyirisquinone, a quinone derivative, can be isolated from the seeds of <i>Iris bungei</i> Maxim [1].
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5604 mL	12.8021 mL	25.6043 mL
5 mM	0.5121 mL	2.5604 mL	5.1209 mL
10 mM	0.256 mL	1.2802 mL	2.5604 mL
50 mM	0.0512 mL	0.256 mL	0.5121 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481