


[Des-Arg10]-HOE 140

Chemical Properties

CAS No. :	1054763-33-5
Formula:	C ₅₃ H ₇₇ N ₁₅ O ₁₂ S
Molecular Weight:	1148.34
Storage:	Keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	[Des-Arg10]-HOE 140 is a bioactive peptide that functions as a potent B1 bradykinin receptor antagonist.
Targets(IC50)	Bradykinin Receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	0.8708 mL	4.3541 mL	8.7082 mL
5 mM	0.1742 mL	0.8708 mL	1.7416 mL
10 mM	0.0871 mL	0.4354 mL	0.8708 mL
50 mM	0.0174 mL	0.0871 mL	0.1742 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481