


(rel)-Tivantinib

Chemical Properties

CAS No. :	905853-99-8
Formula:	C ₂₃ H ₁₉ N ₃ O ₂
Molecular Weight:	369.42
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	(rel)-Tivantinib is a potent, highly selective inhibitor of the receptor tyrosine kinase c-MET and has additional novel targets, GSK3 α and GSK3 β , that are critical in the cellular mechanisms of non-small cell lung cancer (NSCLC) [1].
Targets(IC50)	Others,c-Met/HGFR

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.7069 mL	13.5347 mL	27.0695 mL
5 mM	0.5414 mL	2.7069 mL	5.4139 mL
10 mM	0.2707 mL	1.3535 mL	2.7069 mL
50 mM	0.0541 mL	0.2707 mL	0.5414 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481