

(R)-ICMT-IN-3

Chemical Properties

CAS No. :

Formula:

Molecular Weight:

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	(R)-ICMT-IN-3 (compound ent 2-27) is a potent ICMT inhibitor, with an IC50 value of 0.01 μ M [1].
Targets(IC50)	Others,Transferase

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481