

(R)-2-Isopropyl-8,8,8a-trimethyl-6,7,8,8a-tetrahydrophenanthrene-3,4-diol**Chemical Properties**


CAS No. : 2244560-40-3

Formula: C₂₀H₂₆O₂

Molecular Weight: 298.42

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

**Biological Description**

Description	(R)-2-Isopropyl-8,8,8a-trimethyl-6,7,8,8a-tetrahydrophenanthrene-3,4-diol (compound 48), a catechol compound, is involved in the bioinspired synthesis of rearranged abietane diterpenes associated with pygmaeocins [1].
-------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.351 mL	16.7549 mL	33.5098 mL
5 mM	0.6702 mL	3.351 mL	6.702 mL
10 mM	0.3351 mL	1.6755 mL	3.351 mL
50 mM	0.067 mL	0.3351 mL	0.6702 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286

E_mail: info@targetmol.com

Address: 34 Washington Street, Wellesley Hills, MA 02481