

(Gly22)-amyloid beta-protein(1-40)

Chemical Properties

CAS No. : 175010-18-1

Formula:

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	Amyloid beta-protein(1-40) Arctic variant (Ab40ARC (E22G)), featuring a glycine at position 22 (Gly22), is a peptide utilized in Alzheimer's disease research [1].
Targets(IC50)	Beta Amyloid

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481