

## (Gln22,Asn23)-Amyloid $\beta$ -Protein (1-40)

### Chemical Properties

CAS No. : 374796-75-5

Formula:

Molecular Weight:

Keep away from moisture

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

### Biological Description

Description	(Gln22,Asn23)-Amyloid $\beta$ -Protein (1-40) is a peptide used to study Cerebral Amyloid Angiopathy Mutations [1].
Targets(IC50)	Beta Amyloid

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481