


(E)-p-Coumaryl alcohol

Chemical Properties

CAS No. :	20649-40-5
Formula:	C ₉ H ₁₀ O ₂
Molecular Weight:	150.17
Storage:	Store at low temperature Powder: -20°C for 3 years In solvent: -80°C for 1 year <small>Actual storage temperature shall be subject to the COA.</small>



Biological Description

Description	(E)-p-Coumaryl alcohol ((E)-p-Hydroxycinnamyl alcohol) is a coumarin metabolite isolated from galangal rhizome extract.
Targets(IC50)	Others

Solubility Information

Solubility	DMSO: 100 mg/mL (665.91 mM), Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	6.6591 mL	33.2956 mL	66.5912 mL
5 mM	1.3318 mL	6.6591 mL	13.3182 mL
10 mM	0.6659 mL	3.3296 mL	6.6591 mL
50 mM	0.1332 mL	0.6659 mL	1.3318 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Sukhirun N, et al. Bioefficacy of *Alpinia galanga* (Zingiberaceae) rhizome extracts, (E)-p-acetoxycinnamyl alcohol, and (E)-p-coumaryl alcohol ethyl ether against *Bactrocera dorsalis* (Diptera: Tephritidae) and the impact on detoxification enzyme activities. *J Econ Entomol.* 2011 Oct;104(5):1534-40.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481