

**(±)-Acuminatin****Chemical Properties**


CAS No. : 72881-08-4

Formula: C<sub>21</sub>H<sub>24</sub>O<sub>4</sub>

Molecular Weight: 340.41

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

**Biological Description**

Description	Acuminatin, a novel bis-phenylpropide, has been isolated from Magnolia acuminata L. [1].
-------------	--

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.9376 mL	14.6882 mL	29.3763 mL
5 mM	0.5875 mL	2.9376 mL	5.8753 mL
10 mM	0.2938 mL	1.4688 mL	2.9376 mL
50 mM	0.0588 mL	0.2938 mL	0.5875 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286

E\_mail:info@targetmol.com

Address:34 Washington Street,Wellesley Hills,MA 02481