


**(+)-trans-3'-Acetyl-4'-isobutyrylkhellactone****Chemical Properties**

CAS No. : 1562996-11-5  
 Formula: C<sub>20</sub>H<sub>22</sub>O<sub>7</sub>  
 Molecular Weight: 374.38  
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year  
*Actual storage temperature shall be subject to the COA.*

**Biological Description**

Description	(+)-trans-3'-Acetyl-4'-isobutyrylkhellactone, an angular-type pyranocoumarin, is a natural compound isolated from Peucedani Radix. [1] [2].
-------------	---

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.6711 mL	13.3554 mL	26.7108 mL
5 mM	0.5342 mL	2.6711 mL	5.3422 mL
10 mM	0.2671 mL	1.3355 mL	2.6711 mL
50 mM	0.0534 mL	0.2671 mL	0.5342 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481