


(+)-Saxalin

Chemical Properties

CAS No. :	125517-45-5
Formula:	C ₁₆ H ₁₅ ClO ₅
Molecular Weight:	322.74
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	(+)-Saxalin, a furanocoumarin, is present in Harbouria trachypleura [1].
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.0985 mL	15.4923 mL	30.9847 mL
5 mM	0.6197 mL	3.0985 mL	6.1969 mL
10 mM	0.3098 mL	1.5492 mL	3.0985 mL
50 mM	0.062 mL	0.3098 mL	0.6197 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481