


(+)-Diasyringaresinol**Chemical Properties**

CAS No. : 551-29-1
 Formula: C₂₂H₂₆O₈
 Molecular Weight: 418.44
 Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year
Actual storage temperature shall be subject to the COA.

**Biological Description**

Description	(+)-Diasyringaresinol, a lignan derived from Aquilaria agallocha seeds [1], exhibits potential pharmacological properties.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3898 mL	11.9491 mL	23.8983 mL
5 mM	0.478 mL	2.3898 mL	4.7797 mL
10 mM	0.239 mL	1.1949 mL	2.3898 mL
50 mM	0.0478 mL	0.239 mL	0.478 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481