


3-O-p-Coumaroyloleanolic acid

Chemical Properties

CAS No. :	151334-06-4
Formula:	C ₃₉ H ₅₄ O ₅
Molecular Weight:	602.84
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	3-O-p-Coumaroyloleanolic acid is a natural product for research related to life sciences. The catalog number is TN2976 and the CAS number is 151334-06-4.
Targets(IC50)	Others

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	1.6588 mL	8.2941 mL	16.5881 mL
5 mM	0.3318 mL	1.6588 mL	3.3176 mL
10 mM	0.1659 mL	0.8294 mL	1.6588 mL
50 mM	0.0332 mL	0.1659 mL	0.3318 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel: 781-999-4286 E_mail: info@targetmol.com Address: 34 Washington Street, Wellesley Hills, MA 02481