


## Itazigrel

## Chemical Properties

CAS No. :	70529-35-0
Formula:	C <sub>18</sub> H <sub>14</sub> F <sub>3</sub> N <sub>2</sub> O <sub>2</sub> S
Molecular Weight:	365.37
Storage:	Powder: -20°C for 3 years   In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Itazigrel is a platelet aggregation inhibitor.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.737 mL	13.6848 mL	27.3695 mL
5 mM	0.5474 mL	2.737 mL	5.4739 mL
10 mM	0.2737 mL	1.3685 mL	2.737 mL
50 mM	0.0547 mL	0.2737 mL	0.5474 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

**Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins**

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481