

CM764

Chemical Properties


CAS No. : 1350296-29-5

Formula: C₂₃H₂₇N₄O₃

Molecular Weight: 426.48

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



Biological Description

Description	CM764 is a potent and selective antagonist of sigma-2 receptors, a novel derivative of SN79.
Targets(IC50)	Others,Sigma receptor

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.3448 mL	11.7239 mL	23.4478 mL
5 mM	0.469 mL	2.3448 mL	4.6896 mL
10 mM	0.2345 mL	1.1724 mL	2.3448 mL
50 mM	0.0469 mL	0.2345 mL	0.469 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481