

Compound N025-0038

## Chemical Properties


CAS No. :

Formula: C<sub>16</sub>H<sub>12</sub>O<sub>7</sub>

Molecular Weight: 316.265

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.



## Biological Description

Description	Compound N025-0038 is a useful organic compound for research related to life sciences and the catalog number is T131311.
Targets(IC50)	Others

## Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.1619 mL	15.8093 mL	31.6186 mL
5 mM	0.6324 mL	3.1619 mL	6.3237 mL
10 mM	0.3162 mL	1.5809 mL	3.1619 mL
50 mM	0.0632 mL	0.3162 mL	0.6324 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E\_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481