


Miridesap

Chemical Properties

CAS No. :	224624-80-0
Formula:	C ₁₆ H ₂₄ N ₂ O ₆
Molecular Weight:	340.37
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.



Biological Description

Description	Miridesap (GSK2315698) is a divalent crosslinker of SAP that causes rapid depletion of circulating SAP via hepatic clearance. It also produces significant reduction of SAP from cerebrospinal fluid and visceral amyloid deposits
Targets(IC50)	Beta Amyloid,Others

Solubility Information

Solubility	DMSO: 75 mg/mL (220.35 mM),Sonication is recommended. (< 1 mg/ml refers to the product slightly soluble or insoluble)
In vivo Formulation	10% DMSO+40% PEG300+5% Tween 80+45% Saline: 2 mg/mL (5.88 mM),Sonication is recommended. <i>Please add the solvents sequentially, clarifying the solution as much as possible before adding the next one. Dissolve by heating and/or sonication if necessary. Working solution is recommended to be prepared and used immediately. The formulation provided above is for reference purposes only. In vivo formulations may vary and should be modified based on specific experimental conditions.</i>

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.938 mL	14.6899 mL	29.3798 mL
5 mM	0.5876 mL	2.938 mL	5.876 mL
10 mM	0.2938 mL	1.469 mL	2.938 mL
50 mM	0.0588 mL	0.2938 mL	0.5876 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Reference

Bodin K , Ellmerich S , Kahan M C , et al. Antibodies to human serum amyloid P component eliminate visceral amyloid deposits[J]. NATURE, 2010, 468(7320):93-97.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only· Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481