

Ophiobolin C

Chemical Properties

CAS No. :	19022-51-6
Formula:	C ₂₅ H ₃₈ O ₃
Molecular Weight:	386.57
Storage:	Powder: -20°C for 3 years In solvent: -80°C for 1 year Actual storage temperature shall be subject to the COA.

Biological Description

Description	inhibitor of human CCR5 binding to HIV-1 gp120
Targets(IC50)	Others,HIV Protease,CCR

Solubility Information

Solubility	DMSO: Soluble, Ethanol: Soluble, (< 1 mg/ml refers to the product slightly soluble or insoluble)
------------	---

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5869 mL	12.9343 mL	25.8685 mL
5 mM	0.5174 mL	2.5869 mL	5.1737 mL
10 mM	0.2587 mL	1.2934 mL	2.5869 mL
50 mM	0.0517 mL	0.2587 mL	0.5174 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481