

PPAR γ agonist 2

Chemical Properties

CAS No. :

Formula: C₂₄H₂₀O₅

Molecular Weight: 388.41

Storage: Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Actual storage temperature shall be subject to the COA.

Biological Description

Description	PPAR γ agonist 2 is a highly effective compound that acts as a partial agonist for PPAR γ . It exhibits significant potential for utilization in metabolic disease research [1].
Targets(IC50)	Others,PPAR

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	2.5746 mL	12.873 mL	25.746 mL
5 mM	0.5149 mL	2.5746 mL	5.1492 mL
10 mM	0.2575 mL	1.2873 mL	2.5746 mL
50 mM	0.0515 mL	0.2575 mL	0.5149 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Note: The dilution table applies only to solid products. For liquid products, please calculate the stock solution based on the stated concentration and/or density.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:34 Washington Street,Wellesley Hills,MA 02481